

REF · WEB165

Gold Marsella

Enzymatic Marseille liquid soap



DESCRIPTION

Controlled-foam enzymatic Marseille liquid soap for washing all types of clothes.

DOSAGE AND USE INSTRUCTIONS

Dosage ranges from 5 to 15 g per kg of dry clothes, depending on the degree of dirt. Washing temperature of 40°C is recommended.

TECHNICAL DATA

Appearance Yellowish liquid.

pH 8.5 +/- 0.5

MAIN COMPONENTS

Contains anionic surfactants < 5 %, non-ionic surfactants 5-15 %, soap < 5 %, enzymes < 5 %, optical brighteners < 5 %, perfumes < 5 %, TERPINEOL, LIMONENE, HEXAMETHYLINDANOPYRAN, LAURYLAMINE DIPROPYLENEDIAMINE, BENZISOTHIAZOLINONE, POTASSIUM SORBATE, METHYLCHLOROISOTHIAZOLINONE, METHYLISOTHIAZOLINONE, Alcohols, C12-14(even numbered), ethoxylated(5-15).

PRESENTATION

20 kg, 4x5 kg.

USE OR HANDLING PRECAUTIONS

DANGER



DO NOT SWALLOW

H318 Causes serious eye damage.

P102 Keep out of reach of children.

P262 Do not get in eyes, on skin, or on clothing.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P337+P313 If eye irritation persists: Get medical attention.

Do not mix with other chemicals.

Consult the Safety Data Sheet (SDS) before using the product.

MANUFACTURER

Detervic S.A.

Polígon Industrial Can Malloles · Sant Llorenç Desmunts, 32
08500 Vic · Barcelona

ES · CA · EN · FR · PT

DOCUMENT CONTROL

Rev. 11

Issued: 02/08/2024

Page 1 / 1

SALES CONTACT

T. +34 938 891 344

info@instaquim.com

instaquim.com