

REF · WEB352

Q3

Protein, glycidic or lipidic stain remover.



DESCRIPTION

Q3 is a specific product for removing stains of protein, carbohydrate or lipid origin, such as blood, egg, various sauces, fats, etc. It can be used on most types of laundry, white and coloured, except wool and leather.

DOSAGE AND USE INSTRUCTIONS

It is applied pure directly on the stains. Let work for about 5 min. and rub with a brush if necessary. Then perform a normal wash. Do not exceed contact times of 15 min. A previous test is recommended to ensure the resistance of the colors of the garment. Do not apply heat until the stains have disappeared. Repeat the process if necessary.

TECHNICAL DATA

Appearance Transparent yellow liquid.

pH 8.3 +/- 0.3

MAIN COMPONENTS

Contains phosphonates < 5 %, anionic surfactants < 5 %, non-ionic surfactants 5-15 %, enzymes < 5 %, POTASSIUM SORBATE, BENZISOTHIAZOLINONE, METHYLISOTHIAZOLINONE. Alcohols, C12-14(even numbered), ethoxylated(5-15), C9-C11-alcohol ethoxylated(6), 4-C10-C13-sec-alkylbenzenesulphonic acids.

PRESENTATION

750 ml.

USE OR HANDLING PRECAUTIONS

DANGER



DO NOT SWALLOW

H315

Causes skin irritation.

H318

Causes serious eye damage.

H317

May cause an allergic skin reaction.

P102

Keep out of reach of children.

P280

Wear protective gloves, clothing and eye protection. In case of inadequate ventilation wear respiratory protection.

P303+P361+P353

IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower.

P305+P351+P338

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P501

Dispose of contents/container to hazardous or special waste collection point.

Do not mix with other chemicals.

Consult the Safety Data Sheet (SDS) before using the product.

MANUFACTURER

Detervic S.A.

Polígon Industrial Can Malloles · Sant Llorenç Desmunts, 32
08500 Vic · Barcelona

ES · CA · EN · FR · PT

DOCUMENT CONTROL

Rev. 2

Issued: 25/07/2024

Page 1 / 1

SALES CONTACT

T. +34 938 891 344

info@instaquim.com

instaquim.com